



18101KD, 18102KD Kitchen Stool **Instruction Manual**

Table of Contents

Product Information.....3

Important Safety
Instructions.....4

Stool Assembly.....5

Adjusting The Height.....10

Troubleshooting.....10

Cleaning And
Maintenance.....11

Warranty.....11


Product Information

FEATURES

- For use in kitchen or similar work areas
- Height adjustable legs
- A padded seat and backrest
- Lightweight, comfortable, durable and sturdy
- Australian Standards Approved

Knocks Down

for easy transport



SPECIFICATIONS



Description	Seat Width	Seat Depth	Min Seat Height	Max Seat Height	Width Between Arms	Backrest Height	Max User Weight
Kitchen Stool 18101KD Champagne	360mm	320mm	550mm	710mm	475mm	400mm	110kg
Kitchen Stool 18102KD Greystone	360mm	320mm	550mm	710mm	475mm	400mm	110kg

Important Safety Instructions

READ ALL INSTRUCTIONS BEFORE USING.

WARNING:

1. This Kitchen Stool is designed for a maximum user weight of 110kg. Do not exceed the maximum user weight.
2. Ensure that the leg adjusting pins are fully engaged before use. (See image B). Failure to observe these warnings may result in serious injury.



Kitchen Stool Assembly

A x1 Padded Seat



B x2 Armrests




STEP 1

Connect both armrests (Image B) by sliding them through the tubing holes located on the underside of the seat (Image A).



STEP 2

Ensure the thread on each armrest lines up with the pre-made hole marks on the tubing attached to the base of the seat.



STEP 3

Screw in the four thumbscrews and do not over tighten.





STEP 4

Slide the backrest canes through the tubing located on the underside of the seat.

Note: If the back canes look like they are out of alignment and will not fit properly into the seat base please see the Troubleshooting section on page 10.




STEP 5

There are two safety pins located at the ends of the backrest canes. Press the pin closest to the end in first (Image A), and push to slide the tube into place. Repeat this step for the second pin until you hear the pin click into place. (See Image B)





! Check that all safety pins are engaged, failure to do so may result in serious injury.



Height Adjustment

STEP 1

To adjust the height of the Kitchen Stool, simply press in the locking pin located on the inside of each leg and slide the inner tubing of the leg to the desired height. Ensure the pin slides into one of the holes located on the inner side of the leg.



Troubleshooting

Rework for out of alignment back canes

If the back canes look like they are out of alignment and will not fit properly into the seat base **DO NOT** attempt to force them into alignment. Follow the procedure below or you might damage the seat backrest, screws or canes.

- Loosen (do not remove) the four screws holding the backrest to the back canes.
- Fit the back canes to the base (see Step 4 and 5 in Kitchen Stool Assembly.)
- Re-tighten the four backrest screws.

Cleaning And Maintenance

The Kitchen Stool can be cleaned with a low strength detergent, warm water and a soft textured cloth.

!

Do not use caustic cleaners such as caustic soda or oven cleaner as this may damage the product.

WARRANTY AGAINST DEFECTS IN GOODS

PRODUCT	WARRANTY PERIOD
K CARE TUBALCO R&R	
All Products (Inc. All Electrics)	12 months
Upholstery	12 months
Shower Chairs, Stools & Toilet Seat Raisers (Zinc Treated products)	36 months structural 5 years anti-corrosion (external tube only)
Shower Chairs, Stools & Toilet Seat Raisers (Aluminium Products)	36 months structural 5 years anti-corrosion
Shower Chairs, Stools & Toilet Seat Raisers (Stainless Steel Products)	36 months structural 5 years anti-corrosion
KERRY OXFORD EME EASYVET	
General Equipment	36 months
KERRY	
Manual Handling Equipment	24 months
Slings	12 months
AIR COMFORT	
Frame Compact Lift Chair	36 months
Frame Day Beds	36 months
Foam	36 months
Air Bags and Gas Struts	12 months
Castors	12 months
Fabric	12 months
MY COMFORT	
MC100 & MC310 Frames	10 years
MC100 & MC310 Electrics	2 years

1. K Care Healthcare Solutions provides consumers with the following warranty in relation to this Product, in addition to complying with the requirements of any relevant legislation, including the *Competition and Consumer Act 2010* (Cth) in Australia and the *Consumer Guarantees Act 1993* in New Zealand (the **Acts**), except where a New Zealand consumer acquires the Product for the purposes of a business. The benefits to the consumer given by this warranty are in addition to other rights and remedies of the consumer under the applicable Act.
2. In this warranty, we have used the following definitions:
 - (a) **Form** means a Warranty Claim Form issued by K Care Healthcare Solutions in respect of Products.
 - (b) **K Care Healthcare Solutions** or **our** means K Care Healthcare Solutions Pty Ltd, ABN 47 159 431 099 of 836A Fifteenth Street, Mildura, VIC, 3500.
 - (c) **Products** means the goods manufactured by K Care Healthcare Solutions (including products manufactured by its contract manufacturers) set out in the table above;
 - (d) **Material** means a material or component used by K Care Healthcare Solutions in the manufacture of the Products;
 - (e) **Retailer** means the authorised dealer of Products from whom the Product was purchased;
 - (f) **Warranty Period** means the period within which a defect must appear as set out in the table above, commencing from the date of purchase of the relevant Product. If the Product or part (as the case may be) is repaired or replaced, there will be no extension to the original warranty period; and
 - (g) **Workmanship** means the handling, assembly and manufacturing processes performed by or on behalf of K Care Healthcare Solutions in order to manufacture the Products.
3. K Care Healthcare Solutions warrants that for the duration of the Warranty Period, all Products will be free of faults arising from defects in Workmanship or Materials, on the terms and conditions set out in this warranty.
4. K Care Healthcare Solutions undertakes that if during the Warranty Period any Product, or any part of a Product, has failed to operate correctly due to faulty Workmanship or defective Material, it will repair or replace the Product or part (as the case may be) free of charge (excluding freight fees if applicable) provided that the following procedure is met:
 - (a) The consumer must contact the Retailer or K Care Healthcare Solutions upon becoming aware of any defect to a Product. The consumer will then be provided with a Form, which must be completed by the consumer and returned to K Care Healthcare Solutions at the address, facsimile number or email address stated below together with satisfactory proof of purchase.
 - (b) K Care Healthcare Solutions will review the completed Form to determine whether there is a defect, and if so K Care Healthcare Solutions agrees to (at its option) repair, replace or supply equivalent goods, or pay the cost of any of those remedies to the consumer.
 - (c) If K Care Healthcare Solutions requests the return of the applicable Product or part, the consumer will be responsible for the collection and freight costs of returning that Product or part to K Care Healthcare Solutions.
 - (d) The consumer acknowledges that it is also responsible for the freight costs to deliver any new Product or replacement part to it.
5. To the extent that the supply of Products is a supply to a consumer within the meaning of the applicable Act, the consumer is entitled to consumer guarantees under the applicable Act which K Care Healthcare Solutions does not exclude, restrict or modify. In all other respects, and to the extent permitted by law, K Care Healthcare Solutions:
 - (a) Limits its liability for any non-excludable condition or warranty to rectifying any defect at its option, as set out in paragraph 4(b); and
 - (b) Excludes liability (whether express or implied) of any nature whatsoever for any consequential loss, damage or injury arising as a result of any fault in the Products.
6. This warranty does not extend to damage to Products which occurs during transit or transportation, or which is caused by any abuse, accident or improper installation, connection, use, adjustment or repair or use of goods otherwise than in accordance with instructions issued by K Care Healthcare Solutions.
7. The warranty on Products is waived if any addition or attachment to the Products do not have K Care

Healthcare Solutions approval or are not sold as K Care Healthcare Solutions products. The Products are designed to perform specific tasks under established test loads and unauthorised attachments may produce stresses for which the design is not appropriate.

8. The following applies to consumers who purchased a relevant Product in Australia:
Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

Contact information

K Care Healthcare Solutions	Telephone	1300 783 783
836A Fifteenth Street	Facsimile	1300 784 784
Mildura VIC 3500	Email	customerservice@kcare.com.au

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

All information is considered to be true and correct at the date of publication.
Changes in circumstances after the time of publication may impact on the accuracy
of the information.

202005_1.0

Improving your life everyday™

At Your Service

**Sales, administration
and general enquiries**

National 1300 783 783

Your Customer Service Team

customerservice@kcare.com.au

kcare.com.au