

# Mobility Zone

## Throne 3 in 1 Rail and Spacers



#### Throne 3 in 1

This Rail enables the user to have support close to the body, where it is most needed. This model has a weight loading of 120kg. The user is able to lower safely onto the toilet and pushing up from the seat is also much easier. The Rail is safer than over frames as it is rigidly fixed to the toilet pan, and also does not need to be removed for other users. The handles can also be folded down in order to facilitate a wheelchair transfer. Additionally, either one or both handles can be removed to suit some requirements.



## Technical specifications

- Two optional settings:
  - 510mm or 640mm between the handles and fold down handles
- Height of handles above seat: 260mm approx
- The back plate internal measurements: 370mm

### **Spacers**

The Throne Toilet Spacers have been specifically designed to raise toilet seats. The spacer is supplied with longer bolts to secure the seat and spacer to the pan.

- The Spacers come in two sizes, 50mm and 80mm
- They can be used with or without Throne rails
- Throne spacers allow use of a normal toilet seat and they eliminate movement, which can lead to slipping off the pan.
- Spacers are easily installed, becoming an extension of the toilet pan.

\*\*\*Please note that brackets are required when installing a Throne Rail and an 80mm Spacer.

