



Mobility Zone

EASY GLIDE CHAIR



The Easy Glide has weight activated castors to provide stability when getting in and out of the chair.

- Self-locking castors (release when seated, lock when standing for stability)
- Lightweight aluminium frame
- Adjustments can be made to seat
- Easily cleaned
- Woodgrain look finish
- Available with optional upholstered arm pads
- Large range of fabrics/leather available
- 120kg SWL

Seat Height	Seat width	Seat Depth	Chair Height	Back Height	Chair Width	Arm Height
480mm	420mm	450mm	830mm	370mm	570mm	690mm