



Mobility Zone

Etac Clean Height Adjustable



The Etac Clean Shower Commode is constructed of a powder-coated steel tubing chassis with polypropylene seat/back, arm and foot support. It is secure, safe and simple to use for both carers and users. The height can easily be adjusted, without the need for tools, to suit the user's needs be it for easier transfers or optimal fit over the toilet. The wide push handle makes manoeuvring easy and because it is upright, it ensures a perfect fit over even wall-mounted toilets. The unique foot support has a curved surface for added comfort. The foot support provides side protection and can be pushed under the seat when not in use. There are also alternative foot and calf supports. The open seat at the back facilitates easy access. A selection of accessories are available to enhance comfort, convenience and ensure safe use.

	Item number	Colour	Material	Maximum user weight (kg)	Opening: depth x width (mm)	Seat depth (mm)	Seat height (mm)	Seat width (mm)	Space below seat (mm)	Weight (kg)	Width (mm)	Width between arm supports (mm)
Clean height-adjustable, 4 lockable wheels	80229208	White	Chassis: Powder coated steel tubing. Seat/Back support: Polypropylene. Arm supports: Polypropylene. Foot support: Polypropylene. Castors: Polyamide/TPE/Stainless steel tubing. Castor height adjustment tubes: Stainless steel.	130	360 x 215	440	475 - 600	480	420 - 545*	14,2	540	435

*With pan holders: -15 mm

Powder coated steel frame: Polypropylene seat/back, arm and foot supports

Wide push handle: Upright for better fit over toilets

Gently curved footrest: Can be pushed out of the way under the chair

Chair seat is open at the back: Easy access for hygiene purposes

Height adjustable without the need for tools: Up to 60cm to suit user's needs

Range of positioning, safety, toileting and comfort accessories: 4 Lockable castors

15 Trade Street, Ormiston QLD 4160 Hours: Mon - Fri 8.30 am - 5.00 pm
Sat - 9.00 am - 12.00 pm Phone (07) 3134 4674

www.mobilityzone.com.au