



Rollz Motion 2 - 1 Walker Wheelchair



Features and Benefits:

- Gives you confidence to go further
- Innovative and award-winning design
- Superior brakes suited to wet and dry conditions
- Folds into compact and transportable size
- Strong frame and wheels
- Allows you to stand and walk tall, instead of leaning forward
- Safe with superior manoeuvrability
- Adjustable footrests and handle heights
- Ir-filled tyres to reduce vibration
- Extra padding on the arm rests and handles
- Contoured seat with extra padding.

CODE	DESCRIPTION	UNIT WEIGHT	USER HEIGHT	SEAT HEIGHT	HANDLE HEIGHT	SEAT SIZE WALKER	SEAT SIZE WHEEL CHAIR	OVERALL LENGTH	OVERALL WIDTH	SWL
11101	Rollz Motion - SMALL	11.6kg	1.55m-1.85m	50cm	82-94cm	20cm x 45cm	42cm x 45cm	63cm	67cm	125kg
11102	Rollz Motion - REGULAR	11.6kg	1.6m-1.9m	55cm	87-98cm	20cm x 45cm	42cm x 45cm	63cm	67cm	125kg

15 Trade Street, Ormiston QLD 4160

Mon-Fri 8.30 am - 4.30 pm Sat 9.00 am - 12.00 pm

www.mobilityzone.com.au

07 3134 4674