



Mobility Zone

Jay Fusion Cushion with Fluid Insert



The Jay Fusion Cushion is available in different sizes and with different air inserts.

The Jay Fusion is an advanced new adjustable cushion designed for skin protection. Depending on the user's needs and preferences, this cushion may be ordered with either JAY Flow fluid or ROHO Dry Floatation air inserts, which both perfectly adjust to different sitting positions and ensure proper fluid placement beneath bony prominences to help protect the skin from breakdown. The air in Jay Cushion with Roho Dry Floatation air inserts can be adjusted to accommodate asymmetrical postures and changing needs. The fluid cushion's minimum maintenance design guarantees a long lifespan of the cushion.

The Jay Fusion has a unique dual cover system – the outer cover is made of anti-microbial stretchy machine-washable fabric which reduces surface tension and the inner cover is made of easy to clean, water-resistant fabric.

The cushion provides excellent support for the user's pelvic and thigh area, offers lateral stability for proper positioning and optimal immersion and weight distribution

Key Features

- Adjustable with fluid or air insert options
- Minimum maintenance design
- Incontinent resistant inner cover with moisture resistant zipper
- X-static® Anti-Microbial silver thread stretch cover
- Easy to clean, water resistant fabric

Dimensions

- Width: 14" – 24" (35.5cm – 61cm)
- Depth: 14" - 24" (35.5 – 61cm)
- Height: 4"
- Cushion Weight: 16" x 18" (with fluid insert): 2.4kg
- Cushion Weight: 16" x 18" (with air insert): 1.55kg

Safety

- Maximum user weight: 14" – 21" : 136kg
- Maximum user weight 22" – 24" : 225kg

*More sizes available please contact us for more information

15 Trade Street, Ormiston QLD 4160 Hours: Monday to Friday 8.30am - 5.00pm Phone (07) 3134 4674

www.mobilityzone.com.au

Please contact Aidacare for enquiries relating to sizes outside of the displayed range.